

UNCLASSIFIED

# Update on the Rapid Sample Transfer System for the National Criticality Experiments Research Center

**C. Lu, J. Bounds, E. Sorensen**

LA-UR-14-21789

# Outline

---

- **Introduction**
- **Motivation**
- **Key Requirements**
- **Design**
- **Component Selection**
- **Control**
- **Next Steps**
- **References and Notes**

# Introduction

---

- **Programmatic need to perform measurements on short-lived isotopes in irradiated samples**
- **Samples irradiated at the National Criticality Experiments Research Center (NCERC) in the Device Assembly Facility (DAF) at the Nevada National Security Site (NNSS)**
  - Comet, Flattop, Godiva, Planet
- **Ingress/egress times do not allow for measurement of short-lived isotopes**
- **Pneumatic sample transfer system (rabbit system) is being designed to meet programmatic needs**
  - Rapidly transfer and position samples to and from critical assembly machines and detectors

# Motivation

---

## ■ Scientific:

- Measure short-lived activation products that have not been measured before
- Half lives from seconds to a few hours

## ■ Nuclear Forensics:

- Develop major gamma ray spectrum signatures for rapid post-detonation analysis
- Signatures provide data on intensity of gamma ray emission
- Samples signatures collected post-detonation can be compared with laboratory data to help determine fuel and estimate isotopic ratios and composition

# Motivation

- Net count rates for the top 15 post-fission signatures 4 hours following irradiation:

Natural Uranium

Isotope	Energy (keV)	Net CPS
<sup>134</sup> I	846.77	94.82
<sup>134</sup> I	883.81	62.68
<sup>92</sup> Sr	1383.66	37.34
<sup>142</sup> La	641.07	34.16
<sup>97</sup> Zr	743.03	32.85
<sup>133</sup> I	529.82	31.78
<sup>97</sup> Nb	657.91	29.91
<sup>135</sup> Xe	249.76	28.73
<sup>239</sup> Np	277.766	26.82
<sup>105</sup> Ru	723.97	24.71
<sup>132</sup> Te/ <sup>239</sup> Np	228.22	23.27
<sup>91</sup> Sr	555.47	22.41
<sup>139</sup> Ba	165.96	19.85
<sup>134</sup> Te	210.76	18.62
<sup>105</sup> Ru	676.63	18.40

Highly Enriched Uranium

Isotope	Energy (keV)	Net CPS
<sup>134</sup> I	847.49	172.44
<sup>134</sup> I	884.56	112.75
<sup>92</sup> Sr	1384.2	94.51
<sup>142</sup> La	641.73	74.81
<sup>97</sup> Zr	743.76	68.41
<sup>91</sup> Sr	556.1	58.62
<sup>97</sup> Nb	658.6	58.47
<sup>133</sup> I	530.4	58.28
<sup>135</sup> Xe	250.23	52.99
Annih.	511.48	38.95
<sup>139</sup> Ba	166.27	36.19
<sup>87</sup> Kr	403.14	35.46
<sup>135</sup> I	1260.75	28.24
<sup>88</sup> Kr	196.73	27.96
<sup>85m</sup> Kr	151.4	27.27

UNCLASSIFIED

## Key Requirements - Design

---

- **Transfer capsule must hold samples internally**
- **Minimize potential for activation of transfer capsule**
- **Samples must transfer in less than 3 seconds**
- **Samples can transfer between irradiation point and detector repeatedly**
- **Mechanically adjustable to accommodate any distance from critical assemblies and any detector**
- **Tests must be repeatable**
- **Remote operation from the critical assemblies control room**



UNCLASSIFIED


## Key Requirements – Compressed Air

---

- **DAF house air is compressed air, dried and stored at a nominal 115 psig**
- **Maximum allowable working pressure is 150 psig (DAF design limits)**
  - Rabbit system pressure will be limited to 150 psig by pressure relief valves
- **Rabbit system components will be pressure tested and/or manufacturer certified to operate up to 225 psig**
- **System pressure and transfer capsule speed controlled by a regulator between house air connection and accumulators**

# Design


- Size of the transfer capsule will dictate the minimum bend radius of the tubing
- Size of tubing determines the size of the transfer capsule
- Settled on tubing 1.5 in OD x 1.37 in ID and 1.25 in OD x 3 in length transfer capsule



UNCLASSIFIED

## Design – Transfer Capsule



- 4 designs:



- **Materials will be either high purity polyethylene (PE) or reactor grade graphite**
- **Must build prototypes (smooth and dimpled) to test ease and cost of machining, physical flaws in design, potential catch points, and other unexpected defects**

## Design – Minimum Bend Radius


- Minimum bend radius determined by the dimensions (OD and length) of transfer capsule



$$(R + w)^2 + \left(\frac{L}{2}\right)^2 = (R + D_i)^2$$

- Adding a safety factor of 2 → minimum safe bend radius for nominal dimensions of the transfer capsule = 16.1 in.

# Rabbit System Schematic Diagram



# Component Selection

---

## ■ **Tubing:**

- 1.5 in OD x 1.37 in ID tubing chosen
- Swagelok® sanitary flanges for quick connect/disconnect of sections
- Flanges welded to section ends by Swagelok Southwest Co.

## ■ **Valves:**

- Solenoid operated, check, pressure relief

## ■ **Accumulators**

- Large enough to push transfer capsule ~150 ft in ~1 s

## ■ **Pneumatic linear actuators**

- Adjust transfer station positions

## ■ **Pressure transducers**

## ■ **Position sensors**

- Optical or mechanical





## Component Selection – O-Rings

---

- **2 transfer capsule designs employ o-rings to seal capsule when closed**
- **O-ring materials must have a negligible probability for activation**
  - Many gasket materials have high amounts of F and Cl, both of which have a high  $\sigma_c$
- **Three materials met the requisite mechanical properties and did not contain F or Cl:**
  - Ethylene propylene diene monomer (EPDM)
  - Silicone
  - Medical grade, platinum cured silicone
- **Once received, o-rings were only handled with gloves**
  - The only contamination on the o-rings would be from manufacturer or distributor

# O-Ring Irradiation Experiment

- An irradiation experiment was conducted using Comet on 12 August 2013 to determine the activation probability of each material



# O-Ring Irradiation Experiment

Peak (keV)	Net Count Rate Above Background (cps)	Isotope	Reaction
EPDM			
438.7	1.991 ± 0.026	<sup>69</sup> Zn	<sup>68</sup> Zn(n,γ) <sup>69</sup> Zn
511.0	0.536 ± 0.018	Annih.	
629.9	0.0313 ± 0.0080	<sup>72</sup> Ga	<sup>71</sup> Ga(n,γ) <sup>72</sup> Ga
833.9	0.0667 ± 0.0076	<sup>72</sup> Ga	<sup>71</sup> Ga(n,γ) <sup>72</sup> Ga
1115.4	0.123 ± 0.008	<sup>65</sup> Zn	<sup>64</sup> Zn(n,γ) <sup>65</sup> Zn
1368.7	0.947 ± 0.018	<sup>24</sup> Na	<sup>23</sup> Na(n,γ) <sup>24</sup> Na
1461.0	0.0148 ± 0.0133	<sup>40</sup> K	<sup>39</sup> K(n,γ) <sup>40</sup> K
1525.1	0.0200 ± 0.0037	<sup>42</sup> K	<sup>41</sup> K(n,γ) <sup>42</sup> K
1731.6	0.0358 ± 0.0054	<sup>24</sup> Na	Double Escape
2242.9	0.0716 ± 0.0063	<sup>24</sup> Na	Single Escape
2754.0	0.514 ± 0.013	<sup>24</sup> Na	<sup>23</sup> Na(n,γ) <sup>24</sup> Na
Silicone			
511.0	0.149 ± 0.015	Annih.	
1368.7	2.18 ± 0.03	<sup>24</sup> Na	<sup>23</sup> Na(n,γ) <sup>24</sup> Na
1460.9	0.0185 ± 0.0134	<sup>40</sup> K	<sup>39</sup> K(n,γ) <sup>40</sup> K
1525.1	0.0257 ± 0.0045	<sup>42</sup> K	<sup>41</sup> K(n,γ) <sup>42</sup> K
1731.7	0.0936 ± 0.0073	<sup>24</sup> Na	Double Escape
2242.8	0.172 ± 0.009	<sup>24</sup> Na	Single Escape
2754.0	1.14 ± 0.02	<sup>24</sup> Na	<sup>23</sup> Na(n,γ) <sup>24</sup> Na
Silicone (Platinum Cured)			
2754.0	0.0181 ± 0.0032	<sup>24</sup> Na	<sup>23</sup> Na(n,γ) <sup>24</sup> Na

- Inadequately shielded irradiated Al sample holder determined to be in the line-of-sight of the detector (<sup>24</sup>Na)

# System Control

- NI CompactRIO will hold modules required for control
- Operators will be able to control the system using a LabVIEW interface



# LabVIEW Control Panel

**Rabbit Control Console**

Key Present

Power

Unidirectional  Bidirectional

Triggering

Pressure Select (PSI)

Pressure OK

Capsule at Source Position

Destination Position Empty

Ready to Send

**Source**

Sample Cartridge

Irradiation Position

Detector 1

Detector 2

Detector 3

Detector 4

**Destination**

Irradiation Position

Detector 1

Detector 2

Detector 3

Detector 4

Storage

Capsule at Source Position Override: OFF  ON

Manual Stop

Seconds from Trigger to Arrival at Destination

**Green-when-Empty Positions:**

Irradiation Position

Detector 1

Detector 2

Detector 3

Detector 4

Storage

Sample Cartridge

**Green-when-Good Pressures:**

Accumulator Pressures (PSI) :

Irradiation Position Accumulator

Send Station Accumulator

Count Station Accumulator

**Green-when-Open Valves:**

V-1  V-2  V-3  V-4  V-5  V-6

V-7  V-8  V-9  V-10  V-11  V-12

**Green-when-Extended Linear Actuators:**

LA-1  LA-2  LA-3  LA-4

## Next Steps

---

- Fabricate transfer capsule prototypes
  - No need for high purity PE or reactor grade graphite yet
- Procure components
- Conduct testing in NISC machine shop
  - Prove that we can move the transfer capsule 150 ft in ~1 s
  - Minimum safe bend radius
  - Determine the minimum pressure required to operate
  - Repeatability
- Once the concept has been proven, build the system at DAF

## Acknowledgements

---

This work was supported by the DOE Nuclear Criticality Safety Program, funded and managed by the National Nuclear Security Administration for the Department of Energy.

# References

---

1. Los Alamos National Laboratory, J. Bounds, *Rapid Sample Transfer System at the National Criticality Experiments Research Center (NCERC)*, Memorandum, May 2013. (Unclassified Controlled Nuclear Information Document)
2. MKS Instruments, Inc., *Series 31 ISO-KF Fittings*, 2008
3. Swagelok-TM, *Biopharm Fittings For the Bioprocessing and Pharmaceutical Industries*, Jan 2010
4. The University of Texas at Austin Nuclear Engineering Teaching Laboratory, M. Krause, Personal correspondence, Jun 2013
5. The University of Texas at Austin, K. Berns, J. Fuchs, J. Guidry, M. Yoho, *Automatic Transfer System for Neutron Irradiation Test Bed*, May 2012
6. Agilent Technologies, F. Feyerham, J. Wasson, *The Analysis of Trace Contaminants in High Purity Ethylene and Propylene using GC/MS*, Mar 2005
7. National Aeronautics and Space Administration, A. Teverovsky, *Chlorine contamination diffusion in silicones*, Mar 2002.
8. Parker Hannifin, Email correspondence, Sep 2013.
9. J. H. Chang, "Table of Nuclides," *Korea Atomic Energy Research Institute*, 2000. Website: <http://atom.kaeri.re.kr/>.
10. *Nevada National Security Site Device Assembly Facility Documented Safety Analysis*, 2011. (Unclassified Controlled Nuclear Information Document)
11. *Device Assembly Facility Documented Safety Analysis Addendum for Critical Experiments Facility Operations, Revision 1*, 2010. (Unclassified Controlled Nuclear Information Document)
12. Los Alamos National Laboratory, C. Lu and J. Bounds, *Determination of the Activation Probability of O-Rings for Use in the National Criticality Experiments Research Center Pneumatic Transfer System*, Oct 2013, LA-UR-13-28188.